


VIP-Fixative® pre-filled containers are a safe way to transport, store, and retrieve pathology specimens. Designed with clear polypropylene jars along with a dual-threaded closure system that prevents leakage and reseals time-after-time.

- Containers are leak proof and protect against biohazardous or hazardous material release to the environment
- Range from 7mL dermatology containers to 7.5 Liters operating room organ buckets | filled to 50% capacity
- Container openings range from 13mm to 300mm
- 4 year shelf life | containers marked with a lot code and expiration date
- OSHA compliant labeling
- Manufactured using ISO 13485 and FDA quality systems



Featuring a new label refresh!



VIP-Fixative® Pre-Filled Neutral Buffered Formalin

- primary ingredient is formalin
- fully compatible to replace other 10% NBF fixatives without any adjustments

Available Sizes:

1007	23-730-585	7mL	400 / cs
1015	23-730-586	20mL	100 / cs
1030	23-730-587	40mL	98 / cs
1060	23-730-588	60mL	98 / cs
1090	23-730-589	90mL	72 / cs
1120	23-730-590	120mL	50 / cs
1180	23-730-591	180mL	32 / cs
1480	23-730-579	480mL	9 / cs
2005	23-730-580	960mL	5 / cs
6064	23-730-596	1.9L	6 / cs
6128	23-730-597	3.8L	4 / cs
6500	23-730-577	7.5L	1 / ea

VIP-Fixative® Pre-Filled Zinc Buffered Formalin

- chloride, sulfate, and mercury free
- an excellent replacement for mercury based fixatives

Available Sizes:

1030Z	23-730-598	40mL	98 / cs
1060Z	23-730-599	60mL	98 / cs
1120Z	23-730-600	120mL	50 / cs
1180Z	23-730-601	180mL	32 / cs
1480Z	23-730-602	480mL	9 / cs
2005Z	23-730-603	960mL	5 / cs
6064Z	23-730-604	1.9L	6 / cs
6128Z	23-730-583	3.8L	4 / cs
6500Z	23-730-605	7.5L	1 / ea

* Bulk Sizes Also Available in Neutral and Zinc Buffered Formalin



A complete solution for the treatment, disposal and clean up of formalin / aldehyde waste. Designed as a system for maximum effectiveness in treatment and disposal.

Neutrex products are certified by the California Environmental Protection Agency. The seal indicated that the formalin treatment process results in non-hazardous drain safe waste.



NEUTRA-GUARD® STARTER KIT - 4046
23-730-560

NEUTRA-GUARD® STARTER KIT includes:

- 4 pouches of NEUTRALEX®
- 1 NEUTRA-WIPE® container (60 wipes)
- 1 NEUTRA-FORM®
- 10pk of 8"x10" NEUTRA-PADS®
- 1 NEUTRALEX® pH and Aldehyde Test Kit
- 1 Reaction Container



NEUTRALEX® - 4047
23-730-561

Easy-to-use, pre-measured powder neutralizes formalin in 15 minutes with no drain clogging polymer by-products.



NEUTRA-FORM® - 4304
23-730-569

For accidental spills, one container of control agent will absorb and neutralize a large size spill of 10% formalin or 2% glutaraldehyde.



NEUTRA-PADS® - 4295
23-730-567

Encapsulate and neutralize 10% formalin while reducing vapors during grossing, keeping your area clean and dry.



NEUTRA-WIPES® - 4279
23-730-575

Convenient pop-up dispenser contains Neutra-Wipes® absorb 30-40 mL of 10% NBF and may be used anywhere in the lab.



NEUTRALEX® ALDEHYDE TEST KIT - 4050
23-730-564

pH and Aldehyde Test Kit (ATK) is designed to test pH and residual aldehyde in Neutrex® treated waste.



REACTION CHAMBER - 4276
23-730-565

Waste collection and aldehyde treatment container

9.5 L (2.5 gal) capacity

Neutrex®

4047	23-730-561	Pre Measured Pouches	16 / cs
4048		Bulk 25kg pail	1 / cs
4049		Bulk 9kg pail	1 / ea

Neutra-form®

4304	23-730-569	Neutra-form Container	2 / cs
------	------------	-----------------------	--------

Neutra-pads®

4295	23-730-567	8" x 10"	100 / cs
4295-1	23-730-568	8" x 10"	10 / cs
4297	23-730-570	11" x 17"	15 / cs

Neutra-wipes®


4279	23-730-575	240 wipes / roll	4 / cs
4279K1	23-730-576	60 wipes / roll	1 / cs

Neutra-lex® Aldehyde Test Kit


4050	23-730-564	Test Kit	1 kit / ea
4051		Aldehyde test strips	100 / cs
4052		pH test strips	100 / cs

Neutra-form® Reaction Chamber


4276	23-730-565	Reaction Container	1 / ea
------	------------	--------------------	--------




01a. Scigen 7110 at 40x




01b. Competitor 7211 at 40x




02a. Scigen 7110 at 20x



02b. Competitor 7211 at 20x



03a. Scigen 7110 at 50x



03b. Competitor 7211 at 50x

Scigen® Stains have been designed to be used individually or together as a holistic staining system.

Carefully formulated with pharmaceutical grade materials while adhering to the highest quality standard, ISO 13485, we are able to provide our customers with unparalleled repeatability from batch to batch and proven consistent, reliable results.

Reference our Validation Study for additional information.

Hematoxylin

7110	23-730-561	1 Liter	1 / ea
7111	23-730-561	1 Liter	4 / cs

Eosin-Y

7115	23-730-561	1 Liter	1 / ea
------	------------	---------	--------

Bluing Reagent

7121	23-730-561	1 Liter	1 / ea
------	------------	---------	--------

Differentiator

7119	23-730-614	1 Gallon	1 / ea
------	------------	----------	--------

Stain-Away Wipes

4275	23-730-626	60 Wipes	1 / ea
------	------------	----------	--------



Designed specifically to protect you against 99.9% of germs and bacteria. Formulated with beauty grade ingredients, so your hands will feel smooth and hydrated every time.

Cruelty Free
Silicone and Paraben Free
Synthetic Fragrance Free

Hand Sanitizer | Refreshing Gel Formula

4510		2 oz. w/ flip cap	
4520		4 oz. w/ flip cap	
4525		4 oz. w/ pump cap	12 / cs
4530	23-70-627	8 oz. w/ pump cap	16 / cs
4540	23-70-628	16 oz. w/ pump cap	8 / cs
4550		1 L w/ pump cap	4 / cs
4560	23-70-629	1 gal.	4 / cs

Hand Sanitizer | Spray Formula

4523	23-730-630	4 oz. w/ spray top	12 / cs
------	------------	--------------------	---------



Private label opportunity! Allow us to produce, package and ship your customers’s hand sanitizer needs under their custom label. Great for universities, hospitals, and labs alike.